

Tuff Rock

synthetic resin
Die Stone

Available
in
7 colors

Ivory

Blue

Peach

Topaz

Green

Pink

White



Tuff Rock is scientifically formulated with resins to insure perfect dimensional accuracy every time. It's low expansion of .05% makes it ideal for implants and long span bridgework where precision and accuracy are paramount. Available in 7 pastel colors.



TALLADIUM

27360 W. Muirfield Ln.

Valencia, CA 91355

800-221-6449

web: www.talladium.com

email: info@talladium.com