

Galileo™

DIP POT & BRUSHING WAX



Dip Pot & Wax

Galileo™ Dip Pot has a digital readout that displays the temperature settings and stores it in memory. It heats uniformly and maintains the temperature within one degree. The dimensions of the well are 1 1/8" diameter and 7/8" deep. Dip pot wax uniform dips as thin as 0.2mm with low shrinkage. Also has a high chroma and value making it easy to identify thick or thin areas.

Brushing Wax

Galileo™ Brushing Wax is a unique "Ash Free" wax that can be brushed on as thin as .10mm to cover thin areas on copings, attach sprues to runner bars or base formers. It can also be used to contour pontics. The Brushing Wax is a dead wax with a 40°F lower temperature so it will not create stress and strain when brushed on the parent wax.



TALLADIUM

27360 W. Muirfield Ln.

Valencia, CA 91355

800-221-6449

web: www.talladium.com

email: info@talladium.com