

TALLADIUM, INC. SAFETY DATA SHEET 27360 W. Muirfield Lane, Valencia • California 91355 (661) 295-0900 • (800) 221-6449

Section 1 - Product and Company Identification

Date of SDS: 06/01/1995			Review Date: 02/26/2018
Product Identification: Galileo Waxes:			
Galileo Blockout Wax	Galileo Margin Wax	Galileo Red Utility Wax	Galileo Sculpturing Wax
Galileo Brushing Wax	Galileo Modeling Wax	Galileo Runner Bar Sticks	Galileo Sprue Wax
Galileo Diagnostic Wax	Galileo Pontic Sculpturing Wax	Galileo Runner Bars	Galileo Sticky
Galileo Dip Pot Wax	Galileo Pressable Wax	Galileo Saddle Wax	
Product Synonym: Mixture of Natural Waxes			
Product Class/Use: Dental Wax			
Company Name: Talladium, Incorporated			
Company Address: 27360 West Muirfield Lane			
Valencia, California 91355			
United States of America			
General Information Telephone: (661)-295-0900			
Emergency Telephone: Not required by regulations.			

Section 2 - Composition/Information on Ingredients

These items are Non-Hazardous Materials as defined under the US OSHA Hazard Communication Standard (29CFR 1910.1200)

Section 3 - Hazards Identification, Including Emergency Overview

Fumes from molten product may cause eye or respiratory irritation. Contact with molten product may cause severe thermal burns of both the skin and eyes.

Section 4 - First Aid Measures

Eyes: If contact with wax particles occurs, flush eyes with tepid water for at least 15minutes. If irritation develops, seek medical attention. If contact with molten wax occurs, seek medical attention. Skin: If contact with molten wax, cool wax immediately, but do not remove from skin. Seek medical attention. Inhalation: Remove to fresh air. If breathing becomes difficult, seek medical attention. Ingestion: If small amounts are ingested it will pass the system undigested.

Section 5 - Fire Fighting Measures

There is no explosion hazard from these materials under normal circumstances. However, use of water on molten Wax may cause an explosion.

Extinguishing Media: Dry chemical, foam, carbon dioxide. Do not use water on a wax fire.

Section 6 - Accidental Release Measures

Spill or Leak Procedures: Sweep or vacuum material into waste container for disposal. Let molten was solidify then remove and place into waste container for disposal.

Waste Disposal Method: Contains no hazardous substances, dispose in accordance with Local, State and Federal regulations.

Section 7 - Handling and Storage

Handling: Use caution when handling molten wax Storage: Store in a cool area



TALLADIUM, INC. SAFETY DATA SHEET

Section 8 - Exposure Controls & Personal Protection

Respiratory Protection: None normally required. For fumes use a NIOSH approved organic respirator **Eye Protection:** Use NIOSH approved safety glasses with side shields **Protective Equipment:** N/A

Section 9 - Physical & Chemical Properties

Boiling Point: Not determined **Melting Point:** 75 - 225°F (24 - 107°F) **Appearance:** Solid **Odor:** Slight odor **Solubility:** Not soluble in water

Section 10 - Stability & Reactivity Data

Stability: Stable **Hazardous Decomposition Products:** Carbon monoxide, carbon dioxide, water conditions **Incompatibility:** Strong oxidizing agents **Polymerization:** Will not occur

Section 11 - Toxicological Information

None of the components of these items are damaging.

Section 12 - Ecological Information

Non-hazardous to the environment

Section 13 - Disposal Considerations

Waste Disposal Methods: Follow all Federal, State, and Local regulations.

Section 14 - Transport Information

U.S. DOT Information: Non hazardous, not regulated.

Section 15 - Regulatory Information

These items are non-hazardous

Section 16 - Other Information

The information and recommendations above are believed to be accurate and represent the best information currently available to us. Information prepared to meet U.S. OSHA Hazard Communication Standard, 29 CFR 1910.1200. The person receiving the above material should make their own determination as to the suitability for their use. Talladium makes no warranty expressed or implied with respect to this information and assumes no liability resulting from its use.